



- **L411VB** Laboratory Mixing Faucet, Deck Mounted, 6" Rigid Vacuum Breaker Gooseneck
- **L412VB** Same as Above Except with Swing Vacuum Breaker Gooseneck
- **L414VB** Same as Above Except with Rigid/Swing Vacuum Breaker Gooseneck



NOTES:

1. FIXTURE IS CERTIFIED BY CSA INTERNATIONAL TO COMPLY WITH ASME A112.18.1M AND CSA B125.1.
2. VACUUM BREAKER IS CERTIFIED TO COMPLY WITH ASSE 1001.
3. FIXTURE IS FULLY ASSEMBLED (WITH SUPPLY TUBES LOOSE) AND FACTORY TESTED PRIOR TO SHIPMENT.

Measurements may vary  $\pm 1/4$ ".

© 2008 / WaterSaver Faucet Co.