



# WaterSaver Faucet Co.

Innovative Products for Research

701 West Erie Street Phone 312 666 5500  
Chicago, Illinois 60610 Fax 312 666 8597


## L66WS

DECK MOUNTED DOUBLE SERVICE FIXTURE FOR COLD WATER AND GAS, AIR OR VACUUM WITH MOUNTING SHANKS



WATER SAVER RENEWABLE UNIT WITH REPLACEABLE STAINLESS STEEL SEAT

3/8" IPS MOUNTING SHANKS WITH LOCKNUTS AND WASHERS



NOTE:

CSA CERTIFIED UNDER CAN/CSA B125.M89. COMPLIES WITH ANSI/ASME A112.18.1M.

MEASUREMENTS MAY VARY ±1/4".

Drawing Number: \_\_\_\_\_

Revision Number: 061697-KJS