



WaterSaver


INNOVATIVE PRODUCTS FOR RESEARCH

WaterSaver Faucet
701 W Erie St
Chicago, IL 60654

312 666 5500 TELEPHONE
312 666 5501 FACSIMILE
wsflab.com

L691

PURE WATER FAUCET, TIN-LINED BRASS, DECK MOUNTED, 6" RIGID GOOSENECK



NOTES:

1. FIXTURE HAS INTERIOR LINING OF PURE TIN.
2. FIXTURE IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.

MEASUREMENTS MAY VARY ± 1/4" (6mm).

Drawing Number: _____

Revision Number: 050407-SJP